SLOŽENÍ

40 mg / 2 ml Hyaluronát sodný
30 mg / 2 ml L-prolin
40 mg / 2 ml Glycin
60 mg / 2 ml N-acetylglukosamin

POPIS PRODUKTU

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DESCRIPTION OF DEVICE

Hyaluronic acid is distributed throughout all body tissues and is found in high concentrations in tissues such as the vitreous, synovial fluid, umbilical cord and dermis. Hyaluronic acid acts as a tissue lubricant and is thought to play an important role in regulating interactions between adjacent tissues. It also acts as a viscoelastic support that keeps tissues separate. Different hyaluronic acid preparations may have different molecular weights but are thought to have the same chemical structure.

PURPOSE OF USE AND INDICATIONS

Algyl GGP is intended for the treatment of symptoms of joint disorders such as osteoarthritis.

Algyl GGP is indicated for viscoelastic replenishment, cartilage tissue reconstruction or replacement of synovial fluid in human joints.

CONTRAINDICATIONS

Known sensitivity to any of the ingredients contained in product, pre-existing infections of the skin in the region of the intended injection site, known infection of the index joint, known systemic bleeding disorders. Do not use. infected and inflamed joints.

ADVERSE REACTIONS AND COMPLICATIONS

Infection, arthralgia (knee pain), swelling, joint effusion, stiffness, pain in limb, paraesthesia, phlebitis, pruritus, injection site erythema, injection site edema, injection site pain, injection site reaction, bursitis, immune response.

USER PROFILE

Intraarticular injection should be administered by orthopedist or physiotherapist.

PATIENT PROFILE

Algyl GGP is recommended for patients with osteoarthritis and some other joint diseases. As no clinical evidence is available on the use of hyaluronic acid in children, pregnant and lactating women or in patients with inflammatory joint diseases such as rheumatoid arthritis or Bechterew disease. Algyl GGP is not recommended in these cases.

STORAGE CONDITIONS

The products should be stored at below 25 °C and low humidity (40-60% RH) before being shipped to the customer. Excessive storage conditions such as exposure to direct sunlight and / or storage on or near radiators may destroy the product.
Do not freeze the product.

SINGLE-USE DEVICE - DO NOT RE-USE

The devices are intended for single-use and must not be re-used. Re-use or reprocessing (e.g. cleaning and re-sterilization) may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury, illness or death. Furthermore, reuse or reprocessing of single-use devices may create a risk of contamination e.g. due to the transmission of infectious material from one patient to another. This could result in injury or death of the patient or user. Contaminated implantable substances must not be reprocessed. Any product that has been contaminated by blood, tissue, and/or body fluids/ matter should never be used again and should be handled according to hospital protocol.

STERILIZATION and SHELF-LIFE

Sterilized by moist heat.
The shelf-life is 2 year.

WARNINGS AND PRECAUTIONS


- Those precautions normally considered during injection of substances into joints are recommended.
- Only medical professionals trained in accepted injection techniques for delivering agents to intra-articular synovial joint spaces should inject sodium hyaluronate for this application.
- An excess amount of sodium hyaluronate is not to be used and the patient should be monitored closely.
- The space should not be overfilled.
- If pain increases during the injection procedure, the injection should be stopped and the needle withdrawn.
- Patients experiencing abnormal side effects to the administration of product should consult with a physician immediately.

APPLICATION

- Product is administered once in 6 months.
- Depending on the severity of the disease, the effect of intraarticular injection lasts approximately in 6 months. A second dose can be applied when necessary.
- In case of joint effusion, aspiration, rest, ice pack application and/or intraarticular cortisone administration is recommended to reduce effusion. The preparation can be applied after 2-3 days.
- The contents of the product and the outer surface of the prefilled syringe are sterile until the package is opened.
- Remove the air bubble (if any) before injecting.
- Take the syringe out of the package.

- Attach a suitable cannula (18-22G) carefully by rotating it.

-  Attention
-  See instructions for use
-  Do not re-sterile
-  Lot number
-  Production date
-  Conformity european and notified body number
-  Manufacturer
-  Dry keep
-  Do not re-use
-  Do not use if package is damaged
-  Expiration date
-  Keep away from sunlight
-  Temperature limitation
-  Sterilized using moist heat
-  Sterile fluid path
-  Humidity limitation
-  Non-pyrogenic